

BENDABLE ARMBOARDS



The Dale Bendable ArmBoard is a simple cost effective, tapeless, means to position and secure the wrist to prevent unwanted movement when using arterial and peripheral (intravenous) lines. The ArmBoard is covered in a soft comfortable material and is fastened to the arm with adjustable hook and loop straps. The ArmBoard can be custom shaped to any position.

- limits movement of arterial and other peripheral lines
- easily shaped to desired position
- allows hyperextension of wrist to enhance radial arterial exposure
- adjustable closures and soft material assures fit and comfort


Can be used in the

- Operating Room
- Recovery Room
- ICU
- Patient Care Units

MRI Information

The Dale Bendable ArmBoard is safe to use on patients undergoing MRI under specified conditions.

Large	23.0cm x 9.0cm
Medium	13.7cm x 6.4cm
Small	11.4cm x 4.4cm
X-Small	10.8cm x 2.5cm



LATERAL
MEDICAL

Contact: (07) 3139 1207 or
sales@lateralmedical.com