

# BENDABLE ARMBOARDS



- limits movement of arterial and other peripheral lines
- easily shaped to desired position
- allows hyperextension of wrist to enhance radial arterial exposure
- adjustable closures and soft material assures fit and comfort

## Can be used in the

- Operating Room
- Recovery Room
- ICU
- Patient Care Units

## MRI Information

The Dale Bendable ArmBoard is safe to use on patients undergoing MRI under specified conditions.

The Dale Bendable ArmBoard is a simple cost effective, tapeless, means to position and secure the wrist to prevent unwanted movement when using arterial and peripheral (intravenous) lines.

The ArmBoard is covered in a soft comfortable material and is fastened to the arm with adjustable hook and loop straps. The ArmBoard can be custom shaped to any position.

<b>Large</b>	23.0cm x 9.0cm
<b>Medium</b>	13.7cm x 6.4cm
<b>Small</b>	11.4cm x 4.4cm
<b>X-Small</b>	10.8cm x 2.5cm