

ORAL RAE ENDOTRACHEAL TUBE

Genesis Oral RAE Endotracheal Tube is indicated to be used in any surgery to facilitate directing the tubes away from the surgical field above the neck (ophthalmology, ENT or facial surgery), and also for tonsillectomy and adenoidectomy with the mouth gag.



- Posterior facing bevel with curved atraumatic tip ensures first pass success with reduced laryngeal hang up.
- Designed to for use with video laryngoscopes allowing easy passage.
- Reduced tracheal, laryngeal trauma and vocal cord trauma.

DEVICE CODE	DEVICE NAME	SIZE
DORET80	Oral RAE Endotracheal Tube Cuffed with Posterior Facing Bevel	8.0
DORET75	Oral RAE Endotracheal Tube Cuffed with Posterior Facing Bevel	7.5
DORET70	Oral RAE Endotracheal Tube Cuffed with Posterior Facing Bevel	7.0
DORET65	Oral RAE Endotracheal Tube Cuffed with Posterior Facing Bevel	6.5
DORET60	Oral RAE Endotracheal Tube Cuffed with Posterior Facing Bevel	6.0
DORET55	Oral RAE Endotracheal Tube Cuffed with Posterior Facing Bevel	5.5
DORET50	Oral RAE Endotracheal Tube Cuffed with Posterior Facing Bevel	5.0

Units : Minimum Order Quantity 10 and Multiples thereof