



FIBREOPTIC ORAL/NASAL ENDOTRACHEAL TUBE

- ✓ Designed by clinicians for fibreoptic oral/nasal intubation
- ✓ Suitable for direct laryngoscopy surgical procedures
- ✓ Features extra length and a soft curved atraumatic tip
- ✓ Prevents avoidable naris and larynx trauma
- ✓ Allows for easy securement and manipulation out of surgical field
- ✓ Available in PVC and reinforced in adult lengths
- ✓ Proximal pilot line to prevent trauma to nasal structures
- ✓ Genesis Airway Features
 - Posterior facing bevel
 - Soft, curved, atraumatic tip
 - High volume low pressure cuff



FOIETT and FOIRETT



FOIETT compared to standard ETT

SIZE	PVC	Reinforced
5.0	FOIETT50	FOIRETT50
5.5	FOIETT55	FOIRETT55
6.0	FOIETT60	FOIRETT60
6.5	FOIETT65	FOIRETT65
7.0	FOIETT70	FOIRETT70
7.5	FOIETT75	FOIRETT75
8.0	FOIETT80	FOIRETT80